


Flowflex™

SARS-CoV-2

Antigen Rapid Test (Self-Testing)



A rapid test for the detection of SARS-CoV-2 antigens in nasal specimens. For self-testing use.



The Global Leader of Rapid Test with **26** Years Experience.



Fast



Accurate



Easy to Use



Reliable

ACON®

Flowflex SARS-CoV-2 Antigen Rapid Test (Self-Testing)

The SARS-CoV-2 Antigen Rapid Test is a lateral flow test for the qualitative detection of the antigen from SARS-CoV-2 in nasal swab specimens directly from individuals suspected of COVID-19 within the first seven days of the onset of symptoms. The test can also test specimens from individuals without symptoms.

- **Specimen:** Nasal swab specimens
- **Test Time:** Results at 15 min
- **Shelf life:** 24 months
- **Storage temperature:** 2-30°C
- **Accuracy:** 98.8%
- **Sensitivity:** 97.1%
- **Specificity:** 99.5%


Test Procedure

1. Test Preparation

A. Open your test kit:



B. You should have:




2. Specimen Collection

COLLECTION BY AN ADULT CAREGIVER




OR


SELF COLLECTION



3. Sample Preparation



4. Running the Test



Ordering Information

Product Name	Catalog No.	Format	Specimen	Package
Flowflex SARS-CoV-2 Antigen Rapid Test (Self-Testing)	L031-11855	Cassette	Nasal swabs	1 Test/Kit 20 Kits/Box
	L031-118A5	Cassette	Nasal swabs	5 Tests/Kit



ACON Biotech (Hangzhou) Co., Ltd.
 No.210 Zhenzhong Road, West Lake District, Hangzhou, P.R.China, 310030
 Tel: +86-571-8118 9706 / Fax: +86-571-8777 5781 / Email: information@aconlab.com.cn
www.aconbio.com